



Issue Date: 21st May 2020

SAFETY ALERT NOTICE.

Drug / Product / Device :


“ CROWN NAME FFP2 Mask”

FFP2 Face Mask Marked -

“CROWN NAME”

EN149:2001+A1:2009

CR6920 FFP2 NR



Identified Issue / Concern / Problem/ Evidence of harm.

It appears that the fabric weld which holds the elastic strap in place can be torn when in normal use causing a situation where the mask cannot be secured to the wearers face.

Actions Required / Risk Treatment strategy

Please isolate this specific mask from use, removing from all PPE packs and replacing with another make of FFP2 Mask.

Please ensure that any loose supplied masks of this brand are also withdrawn.

Issued by: Quality, Safety & Risk.